

PABPN1L Rabbit anti-Human Polyclonal (aa163-194) (PE) Antibody - LS-C284474 - LSBio	
CatalogID:	LS-C284474
Target:	poly(A) binding protein, nuclear 1-like (cytoplasmic) (PABPN1L)
Synonyms:	PABPN1L Antibody, EPABP-2 Antibody, PABPNL1 Antibody, EPABP2 Antibody
Host	PABPN1L antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PABPN1L antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PABPN1L antibody was raised against kLH-conjugated synthetic peptide (aa163-194) from internal region of human PABPN1L.
Specificity:	Human PABPN1L
Epitope:	aa163-194
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/30/2014	
© 2014 LifeSpan BioSciences	