

PAX6 Rabbit anti-Human Polyclonal (aa360-390) (PE) Antibody - LS-C284491 - LSBio

CatalogID:	LS-C284491
Target:	paired box 6 (PAX6)
Synonyms:	PAX6 Antibody, AN2 Antibody, D11S812E Antibody, MGDA Antibody, Oculorhombin Antibody, Paired box 6 Antibody, Paired box homeotic gene-6 Antibody, AN Antibody, Aniridia type II protein Antibody, WAGR Antibody, Paired box protein Pax-6 Antibody
Host	PAX6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PAX6 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PAX6 antibody was raised against kLH-conjugated synthetic peptide (aa360-390) of human PAX6.
Specificity:	Human PAX6
Epitope:	aa360-390
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences