

PRKACB Rabbit anti-Human Polyclonal (aa291-319) (FITC) Antibody - LS-C284765 - LSBio

CatalogID:	LS-C284765
Target:	protein kinase, cAMP-dependent, catalytic, beta (PRKACB)
Synonyms:	PRKACB Antibody, Cbeta Antibody, PKACB Antibody, PKA C-beta Antibody, PKA Cbeta Antibody
Family / Subfamily:	Protein Kinase / cAMP
Host	PRKACB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PRKACB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PRKACB antibody was raised against kLH-conjugated synthetic peptide (aa291-319) from the C-terminus of Mouse Prkaca/Prkacb.
Specificity:	Human PRKACB
Epitope:	aa291-319
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences