

PTPLB Rabbit anti-Human Polyclonal (aa224-254) (PE) Antibody - LS-C284875 - LSBio

CatalogID:	LS-C284875
Target:	protein tyrosine phosphatase-like (proline instead of catalytic arginine), member b (PTPLB)
Synonyms:	PTPLB Antibody, HACD2 Antibody
Host	PTPLB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PTPLB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PTPLB antibody was raised against KLH-conjugated synthetic peptide (aa224-254) from the C-terminus of human PTPLB.
Specificity:	Human PTPLB
Epitope:	aa224-254
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/30/2014

© 2014 LifeSpan BioSciences