


ROX / MNT Rabbit anti-Human Polyclonal (aa532-582) Antibody - LS-C290287 - LSBio	
<b>CatalogID:</b>	LS-C290287
<b>Target:</b>	MAX network transcriptional repressor (MNT)
<b>Synonyms:</b>	MNT Antibody, BHLHd3 Antibody, Max-binding protein MNT Antibody, MAX binding protein Antibody, Max-interacting protein Antibody, MNT, MAX dimerization protein Antibody, Protein ROX Antibody, Myc antagonist MNT Antibody, MAD6 Antibody, MXD6 Antibody, ROX Antibody
<b>Host</b>	MNT antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Immunogen Species:</b>	ROX / MNT antibody was raised against Human
<b>Specificity:</b>	Region between residue 532 and 582 of human MAX Binding Protein using the numbering given in entry NP_064706.1 (GeneID 4335).
<b>Epitope:</b>	aa532-582
<b>Reactivity:</b>	Human, Mouse, Rat, Chicken, Xenopus, Zebrafish
<b>Purification:</b>	Immunoaffinity purified
<b>Presentation:</b>	Tris-citrate/phosphate buffer, pH 7-8, 0.09% sodium azide.
<b>Recommended Storage:</b>	Store at 2-8°C for up to 1 year.
<b>Uses:</b>	Western blot (1:2000 - 1:10000), Immunoprecipitation (2 - 10 µg/mg lysate) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Manufacturer:</b>	Bethyl Laboratories, Inc.

**Western Blot Image:**



Detection of Human MNT by Western Blot and Immunoprecipitation. Samples: Whole cell lysate from HeLa (15 and 50 ug for WB; 1 mg for IP, 20% of IP loaded), 293T (T; 50 ug) and Jurkat (J; 50 ug) cells. Antibodies: Affinity purified rabbit anti-MNT antibody used for WB at 0.1 ug/ml (A) and 1 ug/ml (B) and used for IP at 6 ug/mg lysate. MNT was also immunoprecipitated by rabbit anti-MNT antibody A303-626A, which recognizes an upstream epitope. Detection: Chemiluminescence with exposure times of 3 minutes (A) and 30 seconds (B).

**Requested From:**

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences