

CD45 Mouse anti-Rat Monoclonal (MRC OX-30) Antibody - LS-C34503 - LSBio

CatalogID:	LS-C34503
Target:	protein tyrosine phosphatase, receptor type, C (PTPRC)
Synonyms:	PTPRC Antibody, B220 Antibody, CD45R Antibody, CD45 Antibody, gp180 Antibody, Leukocyte common antigen Antibody, L-CA Antibody, Leucocyte common antigen Antibody, LY5 Antibody, T200 Antibody, T200 glycoprotein Antibody, T200 leukocyte common antigen Antibody, CD45 antigen Antibody, LCA Antibody
Family / Subfamily:	Protein Phosphatase / R1/R6
Host	PTPRC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MRC OX-30
Immunogen Species:	CD45 antibody was raised against Rat
Immunogen:	CD45 antibody was raised against lymph node glycoproteins and cells.
Specificity:	Anti-rat CD45 (Leukocyte Common Antigen) monoclonal antibody recognizes a monomorphic determinant of the rat leukocyte common antigen. The antigen recognized is a heavily glycosylated membrane glycoprotein of molecular weight 170000 kD on thymocytes but.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	IHC - Frozen, Flow Cytometry (0.5 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	250 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences