

CD45RC Mouse anti-Rat Monoclonal (MRC OX-22) Antibody - LS-C34516 - LSBio

CatalogID:	LS-C34516
Target:	CD45RC
Host	CD45RC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MRC OX-22
Immunogen Species:	CD45RC antibody was raised against Rat
Immunogen:	CD45RC antibody was raised against phytohemagglutinin-activated rat lymph node cells IMMUNOCYTE.
Specificity:	Rat Leukocyte Common Antigen has been shown to exist in different forms on different lymphoid cell types. Rat L-CA on thymocytes gives one main band at 180 kD, T cells show 4 bands at 180, 190, 200 and 220 kD and B cells; one broad band at 240 kD. T he.
Reactivity:	Rat
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	250 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences