

THY1 / CD90 Mouse anti-Mouse Monoclonal (5a-8) Antibody - LS-C34565 - LSBio

CatalogID:	LS-C34565
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	5a-8
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Immunogen:	THY1 / CD90 antibody was raised against cBA/J.
Specificity:	Purified anti-mouse CD90 (Thy 1.2) is a monoclonal antibody that reacts with all T lymphocytes from mouse strains expressing the Thy 1.2 phenotype (e. g. C57BL/6, C3H/He, DBA/2, CBA/J, BALB/c), but does not react with lymphocytes expressing the Thy 1.1 phe.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	250 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences