

THY1 / CD90 Rat anti-Mouse Monoclonal (IBL-1) Antibody - LS-C34566 - LSBio

CatalogID:	LS-C34566
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	IBL-1
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	THY1 / CD90 antibody was raised against mouse spleen cells.
Specificity:	Anti-mouse CD 90 (Thy 1) monoclonal antibody recognizes the Thy 1 antigen which is a GPI-linked glycoprotein member of the Ig-superfamily. Thy 1, also know as CD90, is expressed on murine T cells, thymocytes, neural cells, Kupffer cells and fibroblasts
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.02% sodium azide.
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	Flow Cytometry (1 µg/10E6 cells) (Optimal dilution to be determined by the researcher)
Size:	250 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences