

B6.1 Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C35202 - LSBio

CatalogID:	LS-C35202
Target:	B6.1 (LOC396098)
Synonyms:	LOC396098 Antibody, B6.1 Antibody, Bu-1b Antibody, BU1B Antibody, B6.2 Antibody, BU1A Antibody
Host	LOC396098 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	B6.1 antibody was raised against Chicken
Immunogen:	B6.1 antibody was raised against bursal cells from one-day-old H. B15 (Bu-1a/b) strain chickens.
Specificity:	Recognizes chicken Bu-1 (chB6), Mr 70-75kD. Does not cross-react with turkey cells when analyzed by immunofluorescent staining and flow cytometry.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunoprecipitation. Immunohistochemistry (acetone-fixed, frozen tissue sections). Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences