

B6.1 Mouse anti-Chicken Monoclonal (FITC) Antibody - LS-C35208 - LSBio

CatalogID:	LS-C35208
Target:	B6.1 (LOC396098)
Synonyms:	LOC396098 Antibody, B6.1 Antibody, Bu-1b Antibody, BU1B Antibody, B6.2 Antibody, BU1A Antibody
Host	LOC396098 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	B6.1 antibody was raised against Chicken
Immunogen:	B6.1 antibody was raised against bursal cells from 2- to 5-week-old CHA (Bu-1b) strain chickens.
Specificity:	Recognizes chicken Bu-1b alloantigen, Mr 70kD. Does not react with cells from CHA and H. B15Ab strains by immunofluorescence.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. Immunoprecipitation. Immunohistochemistry (frozen sections). Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences