

**CD3E Mouse anti-Human Monoclonal (FITC) (UCHT1 ) Antibody - LS-C35367 - LSBio**

<b>CatalogID:</b>	LS-C35367
<b>Target:</b>	CD3e molecule, epsilon (CD3-TCR complex)
<b>Synonyms:</b>	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
<b>Host</b>	CD3E antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	UCHT1
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD3E antibody was raised against Human
<b>Specificity:</b>	Recognizes the 20kD epsilon chain of the CD3/T cell antigen receptor complex.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. In vitro activation. Flow Cytometry: 10 ul/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences