

CD3E Mouse anti-Human Monoclonal (PE, Cy5) Antibody - LS-C35371 - LSBio	
CatalogID:	LS-C35371
Target:	CD3e molecule, epsilon (CD3-TCR complex)
Synonyms:	CD3E Antibody, CD3e antigen Antibody, CD3-epsilon Antibody, T-cell receptor epsilon chain Antibody, TCRE Antibody, T-cell receptors epsilon chain Antibody, T3E Antibody, CD3 Epsilon Antibody
Host	CD3E antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	CD3E antibody was raised against Human
Specificity:	Recognizes the 20kD epsilon chain of the CD3/T cell antigen receptor complex.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Suitable for us in Flow Cytometry, Immunohistochemistry,. Immunoprecipitation and in vitro activation. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: Frozen sections.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	