

CD24 Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C35484 - LSBio	
<b>CatalogID:</b>	LS-C35484
<b>Target:</b>	CD24 molecule
<b>Synonyms:</b>	CD24 Antibody, CD24 molecule Antibody, CD24A Antibody, Signal transducer CD24 Antibody, CD24 antigen Antibody
<b>Host</b>	CD24 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CD24 antibody was raised against Mouse
<b>Specificity:</b>	Recognizes mouse CD24; Mr 95kD.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 4 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	