

CD24 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C35486 - LSBio

CatalogID:	LS-C35486
Target:	CD24 molecule
Synonyms:	CD24 Antibody, CD24 molecule Antibody, CD24A Antibody, Signal transducer CD24 Antibody, CD24 antigen Antibody
Host	CD24 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD24 antibody was raised against Mouse
Specificity:	Recognizes mouse CD24; Mr 95kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Flow Cytometry: 0.3 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences