

CD24 Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C35508 - LSBio	
CatalogID:	LS-C35508
Target:	CD24 molecule
Synonyms:	CD24 Antibody, CD24 molecule Antibody, CD24A Antibody, Signal transducer CD24 Antibody, CD24 antigen Antibody
Host	CD24 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2c,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD24 antibody was raised against Mouse
Specificity:	Recognizes mouse CD24.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry and cell sorting. Immunoprecipitation. Flow Cytometry: 2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	