

CR1 / CD35 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C35560 - LSBio

CatalogID:	LS-C35560
Target:	complement component (3b/4b) receptor 1 (Knops blood group) (CR1)
Synonyms:	CR1 Antibody, C3b/C4b receptor Antibody, C3-binding protein Antibody, CD35 Antibody, C3BR Antibody, Knops blood group antigen Antibody, C4BR Antibody, CD35 antigen Antibody, Complement receptor type 1 Antibody, KN Antibody
Family / Subfamily:	Protease / Serine S1
Host	CR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CR1 / CD35 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CR1 / CD35 antibody was raised against human monocytes.
Specificity:	Recognizes a human complement receptor 1 (CD35) molecule of about 220kD.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 5% glycerol, 0.2% BSA, 0.05% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Suitable for use in Flow Cytometry and Immunohistochemistry. Flow Cytometry: 1:50. Use 80ul.
Uses:	IHC, Flow Cytometry (1:50) (Optimal dilution to be determined by the researcher)
Size:	120 tst
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences