

CD45RA Mouse anti-Human Monoclonal Antibody - LS-C35623 - LSBio	
CatalogID:	LS-C35623
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Immunogen Species:	CD45RA antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RA antibody was raised against stimulated human leukocytes.
Specificity:	Recognizes all CD45 molecules containing the A region exon.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.5 mg/ml Gentamicin sulfate.
Usage Summary:	Suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry.
Uses:	IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	