

CD45RO Mouse anti-Human Monoclonal Antibody - LS-C35628 - LSBio	
<b>CatalogID:</b>	LS-C35628
<b>Target:</b>	CD45RO
<b>Host</b>	CD45RO antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Immunogen Species:</b>	CD45RO antibody was raised against Human
<b>Immunogen:</b>	CD45RO antibody was raised against cultured human T cell line CA-1.
<b>Specificity:</b>	Recognizes the CD45RO isoform of about 180kD.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.5 mg/ml Gentamicin sulfate.
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry.
<b>Uses:</b>	IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 10/1/2014 © 2014 LifeSpan BioSciences	