

CD81 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C35762 - LSBio

CatalogID:	LS-C35762
Target:	CD81 molecule
Synonyms:	CD81 Antibody, CD81 antigen Antibody, CD81 molecule Antibody, CVID6 Antibody, S5.7 Antibody, TAPA1 Antibody, TAPA-1 Antibody, Tetraspanin-28 Antibody, TSPAN28 Antibody, Tspan-28 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD81 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD81 antibody was raised against Human
Immunogen:	CD81 antibody was raised against human B cell line from Burkitt's lymphoma.
Specificity:	Reacts with the 24kD tetraspan CD81 molecule.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	Supplied as a liquid in PBS, pH 7.5, 5% Glycerol, 0.2% BSA, 0.04% sodium azide.
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	Suitable for use in Flow Cytometry and Western Blot.
Uses:	Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences