

CD84 Mouse anti-Human Monoclonal (PE) Antibody - LS-C35772 - LSBio

CatalogID:	LS-C35772
Target:	CD84 molecule
Synonyms:	CD84 Antibody, CD84 antigen Antibody, Cell surface antigen MAX.3 Antibody, CD84 molecule Antibody, HCD84 Antibody, Hly9-beta Antibody, LY9B Antibody, Leukocyte antigen CD84 Antibody, SLAMF5 Antibody, MCD84 Antibody, SLAM family member 5 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD84 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD84 antibody was raised against Human
Immunogen:	CD84 antibody was raised against human MLR stimulated lymphocytes.
Specificity:	Recognizes a molecule of about 74kD.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.5, 15% glycerol, 0.2% BSA, 0.04% sodium azide. Labeled, PE.
Recommended Storage:	May be stored at 4°C. For long-term storage, aliquot and store at 4°C. Do not freeze. Aliquots are stable for at least 12 months. Sensitive to light.
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences