

EPHA2 / EPH Receptor A2 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C36251 - LSBio	
<b>CatalogID:</b>	LS-C36251
<b>Target:</b>	EPH receptor A2 (EPHA2)
<b>Synonyms:</b>	EPHA2 Antibody, CTPA Antibody, EPH receptor A2 Antibody, Ephrin type-A receptor 2 Antibody, Epithelial cell kinase Antibody, ECK Antibody, Sek-2 Antibody, ARCC2 Antibody, Soluble EPHA2 variant 1 Antibody, CTPP1 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Ephrin Receptor
<b>Host</b>	EPHA2 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	EPHA2 / EPH Receptor A2 antibody was raised against Mouse
<b>Immunogen:</b>	EPHA2 / EPH Receptor A2 antibody was raised against mouse EphA2.
<b>Specificity:</b>	Recognizes mouse EphA2.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, BSA.
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	