

LY6G / Gr-1 Rat anti-Mouse Monoclonal (PE, Cy5) (RB6-8C5) Antibody - LS-C36565 - LSBio	
<b>CatalogID:</b>	LS-C36565
<b>Target:</b>	lymphocyte antigen 6 complex, locus G (LY6G)
<b>Synonyms:</b>	LY6G Antibody, Ly-6G Antibody, Gr-1 Antibody, Gr1 Antibody, Ly6g Antibody
<b>Host</b>	LY6G antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Clone Name:</b>	RB6-8C5
<b>Conjugations:</b>	Phycoerythrin (PE), Cy5
<b>Immunogen Species:</b>	LY6G / Gr-1 antibody was raised against Mouse
<b>Specificity:</b>	Recognizes mouse Ly-6G, Mr 21-25kD.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Western blotting. In vitro complement-dependent cytotoxicity assays. In vivo depletion of Ly-6G+ cells. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 10/1/2014 © 2014 LifeSpan BioSciences	