

LY6G / Gr-1 Rat anti-Mouse Monoclonal (PE, Cy5.5) (RB6-8C5) Antibody - LS-C36568 - LSBio

CatalogID:	LS-C36568
Target:	lymphocyte antigen 6 complex, locus G (LY6G)
Synonyms:	LY6G Antibody, Ly-6G Antibody, Gr-1 Antibody, Gr1 Antibody, Ly6g Antibody
Host	LY6G antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	RB6-8C5
Conjugations:	Phycoerythrin (PE), Cy5.5
Immunogen Species:	LY6G / Gr-1 antibody was raised against Mouse
Antigen Type:	Protein
Immunogen:	LY6G / Gr-1 antibody was raised against mouse Gr-1.
Specificity:	Recognizes mouse Ly-6G, Mr 21-25kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Suitable for use in Flow Cytometry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Western Blot, In-vitro complement-dependent cytotoxicity assays, and In-vivo depletion of Ly-6G+ cells. Flow Cytometry: 0.2 ug/1E6 cells. Immunohistochemistry: frozen sections.
Uses:	IHC - Frozen, Immunofluorescence, Western blot, Immunoprecipitation, Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences