

LY6G / Gr-1 Rat anti-Mouse Monoclonal (PE, Cy7) (RB6-8C5) Antibody - LS-C36569 - LSBio	
CatalogID:	LS-C36569
Target:	lymphocyte antigen 6 complex, locus G (LY6G)
Synonyms:	LY6G Antibody, Ly-6G Antibody, Gr-1 Antibody, Gr1 Antibody, Ly6g Antibody
Host	LY6G antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	RB6-8C5
Conjugations:	Phycoerythrin (PE), Cy7
Immunogen Species:	LY6G / Gr-1 antibody was raised against Mouse
Specificity:	Recognizes mouse Ly-6G, Mr 21-25kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Western blotting. In vitro complement-dependent cytotoxicity assays. In vivo depletion of Ly-6G+ cells. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	