

HIV-1 env Mouse anti-Human Immunodeficiency Virus Monoclonal (FITC) Antibody - LS-C36722 - LSBio	
CatalogID:	LS-C36722
Target:	gp160; envelope glycoprotein
Synonyms:	env Antibody, HIV-1 Envelope Glycoprotein Antibody
Host	env antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HIV-1 env antibody was raised against Human Immunodeficiency Virus
Specificity:	Binds to native and recombinant HIV-1 gp120 in Dot Blot assays and stains cell surfaces in direct immunofluorescence assays.
Reactivity:	Human Immunodeficiency Virus
Purification:	Purified
Presentation:	10 mM Tris, pH 8.0. No stabilizers or preservatives added.
Usage Summary:	Direct Immunofluorescent staining of cells may be performed in 1:50 dilution range. Dot Blot assays with purified gp120 may be performed in 1:200 dilution range.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/1/2014 © 2014 LifeSpan BioSciences	