

IL17C Mouse anti-Human Monoclonal (CF) Antibody - LS-C37040 - LSBio

CatalogID:	LS-C37040
Target:	interleukin 17C (IL17C)
Synonyms:	IL17C Antibody, CX2 Antibody, Cytokine CX2 Antibody, Interleukin 17C Antibody, IL-17C Antibody, IL-21 Antibody, Interleukin-17C Antibody
Family / Subfamily:	Interleukin
Host	IL17C antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Carboxyfluorescein (CF)
Immunogen Species:	IL17C antibody was raised against Human
Immunogen:	IL17C antibody was raised against human IL-17C.
Specificity:	Recognizes human IL-17C.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS.
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences