

FCRH1 / FCRL1 Mouse anti-Human Monoclonal (PE) Antibody - LS-C37090 - LSBio

CatalogID:	LS-C37090
Target:	Fc receptor-like 1 (FCRL1)
Synonyms:	FCRL1 Antibody, CD307a Antibody, CD307a antigen Antibody, Fc receptor-like 1 Antibody, FcR-like protein 1 Antibody, Fc receptor-like protein 1 Antibody, FCRH1 Antibody, IFGP1 Antibody, IRTA5 Antibody, HIFGP1 Antibody, IFGP family protein 1 Antibody, Fc receptor homolog 1 Antibody
Host	FCRL1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FCRH1 / FCRL1 antibody was raised against Human
Immunogen:	FCRH1 / FCRL1 antibody was raised against human FcRH1/IRTA5.
Specificity:	Recognizes human FcRH1/IRTA5.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences