

Ly6C Rat anti-Mouse Monoclonal (FITC) (HK1.4) Antibody - LS-C37237 - LSBio

CatalogID:	LS-C37237
Target:	Ly6-C antigen
Host	Ly6-C antigen antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2c,k
Clone Name:	HK1.4
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Ly6C antibody was raised against Mouse
Specificity:	Recognizes murine Ly-6C, Mr 18kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. In vitro activation of T lymphocytes. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences