

NKG2A/C/E Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C37669 - LSBio

CatalogID:	LS-C37669
Target:	NKG2A/C/E
Host	NKG2A/C/E antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NKG2A/C/E antibody was raised against Mouse
Immunogen:	NKG2A/C/E antibody was raised against CHO cells transfected with murine NKG2.
Specificity:	Recognizes mouse NKG2 (Mr ~43kD).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences