

Klr1c / CD161 Mouse anti-Mouse Monoclonal (PE) (PK136) Antibody - LS-C37695 - LSBio	
<b>CatalogID:</b>	LS-C37695
<b>Target:</b>	killer cell lectin-like receptor subfamily B member 1C (Klr1c)
<b>Synonyms:</b>	Klr1c Antibody, AI462337 Antibody, CD161 Antibody, KLRBC Antibody, Ly-59 Antibody, Ly59 Antibody, NK-1 Antibody, Nk-1.2 Antibody, NK1 Antibody, Nk1.2 Antibody, NKRP1 Antibody, NK-RP1 Antibody, NK1.1 Antibody, NKRP140 Antibody, Ly-55c Antibody, Ly55c Antibody, NKR-P1.9 Antibody, Nkrp1c Antibody
<b>Host</b>	Klr1c antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	PK136
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	Klr1c / CD161 antibody was raised against Mouse
<b>Immunogen:</b>	Klr1c / CD161 antibody was raised against cD mouse spleen and bone marrow cells enriched for NK-1+ cells.
<b>Specificity:</b>	Recognizes mouse NK1.1 alloantigen, Mr 80-85kD (39kD under reducing conditions).
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Usage Summary:</b>	Flow cytometry. In vitro complement-dependent cytotoxicity assays. Immunoprecipitation. In vivo depletion and in vitro activation of NK1.1 cells. In vitro redirected lysis. Flow Cytometry: 0.5 ug/1E6 cells.
<b>Uses:</b>	Immunoprecipitation, Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 10/1/2014 © 2014 LifeSpan BioSciences	