

DEFA5 / Defensin 5 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C37788 - LSBio	
CatalogID:	LS-C37788
Target:	defensin, alpha 5, Paneth cell-specific (DEFA5)
Synonyms:	DEFA5 Antibody, DEF5 Antibody, Defensin, alpha 5 Antibody, Defensin 5 Antibody, Defensin-5 Antibody, HD-5 Antibody, HD5(20-94) Antibody
Host	DEFA5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DEFA5 / Defensin 5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DEFA5 / Defensin 5 antibody was raised against a synthetic17-aa peptide sequence near the C-terminus of human mature alpha defensin-5 (KLH coupled).
Specificity:	Recognizes human alpha defensin-5. No significant homology exists with other alpha-defensins.
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Suitable for use in ELISA and Western Blot. Western Blot: 1-10 ug/ml using ECL. ELISA: 0.5-1 ug/ml.
Uses:	Western blot (1 - 10 µg/ml), ELISA (0.5 - 1 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	