

SLAMF6 / NTBA Mouse anti-Human Monoclonal (CF) Antibody - LS-C37813 - LSBio

CatalogID:	LS-C37813
Target:	SLAM family member 6 (SLAMF6)
Synonyms:	SLAMF6 Antibody, Activating NK receptor Antibody, CD352 Antibody, CD352 antigen Antibody, NK-T-B-antigen Antibody, NTB-A Antibody, NTBA receptor Antibody, KAL1b Antibody, SF2000 Antibody, SLAM family member 6 Antibody, Ly108 Antibody, NTBA Antibody, KALI Antibody
Host	SLAMF6 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Carboxyfluorescein (CF)
Immunogen Species:	SLAMF6 / NTBA antibody was raised against Human
Immunogen:	SLAMF6 / NTBA antibody was raised against human NTB-A/SLAMF6.
Specificity:	Recognizes human NTB-A/SLAMF6.
Reactivity:	Human
Purification:	Purified
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences