

CXCL4 / PF4 Mouse anti-Human Monoclonal (CF) Antibody - LS-C38158 - LSBio

CatalogID:	LS-C38158
Target:	platelet factor 4 (PF4)
Synonyms:	PF4 Antibody, C-X-C motif chemokine 4 Antibody, CXCL4 Antibody, Iroplact Antibody, Oncostatin-A Antibody, PF-4 Antibody, SCYB4 Antibody, Platelet factor 4 Antibody
Family / Subfamily:	Interocrine / Interocrine alpha
Host	PF4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Carboxyfluorescein (CF)
Immunogen Species:	CXCL4 / PF4 antibody was raised against Human
Immunogen:	CXCL4 / PF4 antibody was raised against human CXCL4/PF4.
Specificity:	Recognizes human CXCL4/PF4.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	Saline, 0.5% BSA, 0.1% sodium azide.
Recommended Storage:	+4°C (short-term), -20°C (long-term). Protect from light.
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: 10 ul label 5x10 ⁵ cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences