

PDGFRA / PDGFR Alpha Mouse anti-Human Monoclonal (CF) Antibody - LS-C38270 - LSBio	
<b>CatalogID:</b>	LS-C38270
<b>Target:</b>	platelet-derived growth factor receptor, alpha polypeptide (PDGFRA)
<b>Synonyms:</b>	PDGFRA Antibody, CD140A Antibody, CD140a antigen Antibody, PDGFalphaR Antibody, PDGFR2 Antibody, PDGFRalpha Antibody, AlphaPDGFR Antibody, RHEPDGFRA Antibody, PDGF receptor alpha Antibody, PDGFR-2 Antibody, PDGFRA/BCR fusion Antibody, APDGFR Antibody, Pdgf alpha-receptor Antibody, PDGF-R-alpha Antibody, PDGFR-alpha Antibody
<b>Family / Subfamily:</b>	Protein Kinase / PDGF Receptor
<b>Host</b>	PDGFRA antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Conjugations:</b>	Carboxyfluorescein (CF)
<b>Immunogen Species:</b>	PDGFRA / PDGFR Alpha antibody was raised against Human
<b>Immunogen:</b>	PDGFRA / PDGFR Alpha antibody was raised against human PDGF Ra.
<b>Specificity:</b>	Recognizes human PDGF Ra.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, BSA.
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	