

PDGFRB / PDGFR Beta Mouse anti-Human Monoclonal (CF) Antibody - LS-C38281 - LSBio	
CatalogID:	LS-C38281
Target:	platelet-derived growth factor receptor, beta polypeptide (PDGFRB)
Synonyms:	PDGFRB Antibody, CD140B Antibody, Beta pdgfr Antibody, Beta type PDGFR Antibody, BPDGFR Antibody, CD140b antigen Antibody, IBGC4 Antibody, JTK12 Antibody, PDGF-beta R Antibody, PDGF-R-beta Antibody, PDGFR Antibody, PDGFRbeta Antibody, PDGF receptor beta Antibody, PDGF-R Antibody, PDGFbetaR Antibody, PDGFR Beta Antibody, PDGFR-1 Antibody, PDGFR-beta Antibody, BetaPDGFR Antibody, Pdgf beta-receptor Antibody, PDGFR1 Antibody, PDGFRB/CEV14 fusion gene Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	PDGFRB antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Carboxyfluorescein (CF)
Immunogen Species:	PDGFRB / PDGFR Beta antibody was raised against Human
Immunogen:	PDGFRB / PDGFR Beta antibody was raised against human PDGF Rb.
Specificity:	Recognizes human PGDF Rb.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.5% BSA, 0.1% sodium azide.
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: 10 ul labels 1x1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	