

STAT3 Mouse anti-Human Monoclonal (CF) Antibody - LS-C39264 - LSBio

CatalogID:	LS-C39264
Target:	signal transducer and activator of transcription 3 (acute-phase response factor) (STAT3)
Synonyms:	STAT3 Antibody, Acute-phase response factor Antibody, APRF Antibody, DNA-binding protein APRF Antibody, HIES Antibody
Host	STAT3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Carboxyfluorescein (CF)
Immunogen Species:	STAT3 antibody was raised against Human
Immunogen:	STAT3 antibody was raised against human STAT3.
Specificity:	Recognizes human Stat3.
Reactivity:	Human
Purification:	Purified
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry:.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences