

CD7 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C3959 - LSBio

CatalogID:	LS-C3959
Target:	CD7 molecule
Synonyms:	CD7 Antibody, CD7 antigen Antibody, CD7 antigen (p41) Antibody, CD7 molecule Antibody, LEU-9 Antibody, T-cell leukemia antigen Antibody, T-cell antigen CD7 Antibody, T-cell surface antigen Leu-9 Antibody, TP41 Antibody, gp40 Antibody, p41 protein Antibody, Tp40 Antibody
Host	CD7 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD7 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD7 antibody was raised against normal human peripheral blood lymphocytes.
Specificity:	CD7 is a single chain cell surface glycoprotein of 40kD, expressed by T cells, NK cells and thymocytes. CD7 is also expressed by tumor cells in some cases of acute myeloid leukaemia.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.1% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Flow Cytometry: Neat.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences