

**CD8B / CD8 Beta Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C4063 - LSBio**

<b>CatalogID:</b>	LS-C4063
<b>Target:</b>	CD8b molecule
<b>Synonyms:</b>	CD8B Antibody, CD8b antigen Antibody, CD8b molecule Antibody, CD8B1 Antibody, p37 Antibody, LY3 Antibody, LYT3 Antibody, Leu2 Antibody
<b>Host</b>	CD8B antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD8B / CD8 Beta antibody was raised against Chicken
<b>Specificity:</b>	Recognizes the beta chain (Mr 34kD) of chicken CD8.
<b>Reactivity:</b>	Chicken
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (acetone-fixed, frozen tissue sections). Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences