

CD9 Mouse anti-Human Monoclonal (PE) Antibody - LS-C4078 - LSBio	
CatalogID:	LS-C4078
Target:	CD9 molecule
Synonyms:	CD9 Antibody, 5H9 antigen Antibody, BA2 Antibody, BA-2/p24 antigen Antibody, CD9 antigen Antibody, CD9 antigen (p24) Antibody, DRAP-27 Antibody, MIC3 Antibody, MRP-1 Antibody, Leukocyte antigen MIC3 Antibody, Tetraspanin-29 Antibody, TSPAN29 Antibody, BTCC-1 Antibody, CD9 molecule Antibody, Motility related protein-1 Antibody, Motility-related protein Antibody, TSPAN-29 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD9 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD9 antibody was raised against Human
Specificity:	Recognizes human CD9.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Western blotting. Platelet inhibition and aggregation studies. Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Western blot, Immunoprecipitation, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	