


BrdU Mouse Monoclonal (MoBU-1) Antibody - LS-C40882 - LSBio

CatalogID:	LS-C40882
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MoBU-1
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	+4°C, avoid freezing
Uses:	IHC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Uses Not Recommended:	IHC - Paraffin
Size:	100 µg
Concentration:	0.5 mg/ml

Flow Cytometry Image:



BrdU-incorporated Hut-78 cells stained with purified MoBU-1, followed by anti-mouse IgG FITC

Requested From:	Japan
------------------------	-------

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences