


**CD45RO Mouse anti-Human Monoclonal (UCHL1) Antibody - LS-C41089 - LSBio**

<b>CatalogID:</b>	LS-C41089
<b>Target:</b>	CD45RO
<b>Host</b>	CD45RO antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	UCHL1
<b>Immunogen Species:</b>	CD45RO antibody was raised against Human
<b>Reactivity:</b>	Human, Chimpanzee
<b>Purification:</b>	Affinity purified
<b>Recommended Storage:</b>	+4°C, avoid freezing
<b>Uses:</b>	IHC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg

**Flow Cytometry Image:**



Human peripheral blood lymphocytes stained with purified UCHL1, followed by anti-mouse IgGs FITC

**Requested From:** Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences