

**Complement C1QA Sheep anti-Human Polyclonal (FITC) Antibody - LS-C42308 - LSBio**

<b>CatalogID:</b>	LS-C42308
<b>Target:</b>	complement component 1, q subcomponent, A chain (C1QA)
<b>Synonyms:</b>	C1QA Antibody, Complement C1QA Antibody
<b>Host</b>	C1QA antibody was produced in Sheep
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Complement C1QA antibody was raised against Human
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Complement C1QA antibody was raised against human C1q, purified from plasma.
<b>Specificity:</b>	Recognizes human C1q, a sub-component of the C1 complement component
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Uses:</b>	IHC - Frozen (Optimal dilution to be determined by the researcher)
<b>Size:</b>	1,000 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences