

DPP4 / CD26 Mouse anti-Human Monoclonal (RPE) Antibody - LS-C43463 - LSBio

CatalogID:	LS-C43463
Target:	dipeptidyl-peptidase 4 (DPP4)
Synonyms:	DPP4 Antibody, ADCP-2 Antibody, ADABP Antibody, Dipeptidyl peptidase IV Antibody, Dipeptidyl-peptidase 4 Antibody, DPPIV Antibody, ADCP2 Antibody, T-cell activation antigen CD26 Antibody, CD26 Antibody, CD26 antigen Antibody, Dipeptidyl peptidase 4 Antibody, Dipeptidylpeptidase 4 Antibody, DPP IV Antibody, TP103 Antibody
Family / Subfamily:	Protease / Serine S9B
Host	DPP4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	DPP4 / CD26 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	DPP4 / CD26 antibody was raised against human T-cell leukaemia cells (T-CLL).
Specificity:	Recognizes the human CD26 cell surface antigen, a 110kD molecule expressed by activated T cells, B cells and macrophages. This molecule is also known as dipeptidylpeptidase IV
Reactivity:	Human
Purification:	Purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, PBS, pH 7.4, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells or 100 ul whole blood. Method sheets are available on request.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences