

TNF Alpha Rat anti-Human Monoclonal (RPE) Antibody - LS-C43750 - LSBio	
CatalogID:	LS-C43750
Target:	tumor necrosis factor (TNF)
Synonyms:	TNF Antibody, APC1 protein Antibody, Cachectin Antibody, DIF Antibody, TNF Alpha Antibody, Tumor necrosis factor-alpha Antibody, TNF-alpha Antibody, TNF-a Antibody, TNFA Antibody, Tumor necrosis factor Antibody, TNF, macrophage-derived Antibody, TNF, monocyte-derived Antibody, TNFSF2 Antibody
Family / Subfamily:	TNF / not assigned-TNF
Host	TNF antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	TNF Alpha antibody was raised against Human
Reactivity:	Human
Purification:	Purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, PBS, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 1×10^6 cells in 100 ul. 9alphavers=NA">BUF09<a/>) for this purpose.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/1/2014	
© 2014 LifeSpan BioSciences	