

CD55 Mouse anti-Human Monoclonal (RPE) Antibody - LS-C43794 - LSBio

CatalogID:	LS-C43794
Target:	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
Synonyms:	CD55 Antibody, CD55 antigen Antibody, Decay-accelerating factor Antibody, TC Antibody, CR Antibody, CROM Antibody, DAF Antibody
Host	CD55 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD55 antibody was raised against Human
Reactivity:	Human
Purification:	Protein G purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, PBS, 0.09% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences