

CD44 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C43853 - LSBio

CatalogID:	LS-C43853
Target:	CD44 molecule (Indian blood group)
Synonyms:	CD44 Antibody, CDW44 Antibody, CD44 antigen Antibody, Cell surface glycoprotein CD44 Antibody, CSPG8 Antibody, ECMR-III Antibody, Epican Antibody, HUTCH-I Antibody, LHR Antibody, MDU2 Antibody, MC56 Antibody, MIC4 Antibody, Hermes antigen Antibody, Hyaluronate receptor Antibody, Phagocytic glycoprotein 1 Antibody, PGP-1 Antibody, PGP-I Antibody, Pgp1 Antibody, Phagocytic glycoprotein I Antibody, MDU3 Antibody, CD44R Antibody, HCELL Antibody, Heparan sulfate proteoglycan Antibody
Host	CD44 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD44 antibody was raised against Human
Reactivity:	Human
Purification:	Purified
Presentation:	Tris buffered saline, 0.02% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences