



CD4 Antigen (CD4) Mouse Monoclonal (RPE) Antibody

<b>CatalogID:</b>	LS-C43875
<b>Target Protein:</b>	CD4 Antigen (CD4)
<b>Gene Synonyms:</b>	CD4, CD4MUT, t-cell surface glycoprotein cd4
<b>Host:</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b
<b>Modifications:</b>	R. Phycoerythrin (RPE)
<b>Antigen Species:</b>	
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	Cells
<b>Specificity:</b>	Leucocytes isolated from porcine gut lamina propria
<b>Species Reactivity:</b>	Pig
<b>Purification:</b>	Protein G column
<b>Volume:</b>	100 $\mu$ l
<b>Presentation:</b>	Phosphate buffered saline pH7.4
<b>Uses:</b>	IHC-Fr (Optimal dilution to be determined by the researcher)
<b>Summary:</b>	Recognises the porcine CD4a cell surface antigen.

Laboratory Reagent For In Vitro Research Use Only

Created on 8/21/2008

© 2008 LifeSpan BioSciences