

CD45RA Mouse anti-Pig Monoclonal (RPE) (MIL13) Antibody - LS-C43879 - LSBio

CatalogID:	LS-C43879
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MIL13
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD45RA antibody was raised against Pig
Immunogen:	CD45RA antibody was raised against cells isolated from porcine mesenteric lymph node.
Specificity:	Recognizes the porcine CD45RA cell surface antigen
Reactivity:	Pig
Purification:	Purified
Reconstitution:	Distilled Water.
Presentation:	Lyophilized, 0.09% sodium azide, 1% BSA.
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences