

CD45RA Mouse anti-Dog Monoclonal Antibody - LS-C43971 - LSBio

CatalogID:	LS-C43971
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	CD45RA antibody was raised against Dog
Reactivity:	Dog
Presentation:	Tissue culture supernatant, 0.09% sodium azide.
Usage Summary:	Flow Cytometry: Use 10 ul of the suggested working dilution to label 10 ⁶ cells in 100 ul. Method sheets are available on request.
Uses:	IHC - Frozen, Flow Cytometry (1:1) (Optimal dilution to be determined by the researcher)
Size:	2,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences